

★  
Excellence by  
**Contour Group**  
contourgroup.co.uk



# CONTOUR DERMA II

✦ ✦ *REVITALISE & REJUVENATE YOUR CLIENT'S SKIN* ✦ ✦

NON-SURGICAL ADVANCED SKIN REMODELING & TIGHTENING DEVICE



**POTENTIAL  
EARNINGS:  
£300-£600  
PER HOUR!**

# CONTOUR DERMA II

## RADIO FREQUENCY + MICRO-NEEDLING

### TARGET & REPAIR AGEING SKIN ISSUES ON THE FACE & BODY SKIN REJUVENATION // FRACTIONAL REMODELING



#### TECHNICAL SPECIFICATIONS

RF Frequency: 1MHZ

Top Energy: 35W

Tips: Invasive

Tip Size: 15\*15mm

Tip Pin: 25 & 64 pins

Needle Size: Diameter 0.12mm

Needle Pin Depth: 0.5mm - 3.5mm

Screen Size: 10.4 inch touch screen

Dimensions: 50cm x 50cm x 110cm

Net Weight: 40kg



ULTIMATE SKIN RENEWAL



TIGHTENING & LIFTING



COLLAGEN & ELASTIN INDUCTION



IMPROVES SKIN LAXITY, TONE & TEXTURE



FINE LINE & WRINKLE REDUCTION



TARGETS COMMON SKIN CONDITIONS



TREATS STRETCH MARKS, SCARS & BURNED TISSUE

Utilise this **advanced skin rejuvenation system** which causes mild induced trauma stimulating the natural healing process to help clients achieve a reduced appearance of stretch marks, wrinkles, and fine lines, as well as reducing pore size and improving skin tone and texture. Additionally, reduce pigmentation to attain glowing, radiant, and renewed skin.


★ **SCIENTIFICALLY PROVEN TO PRODUCE FASTER & MORE  
DRAMATIC RESULTS OVER TRADITIONAL MICRO-NEEDLING**

# RENEW & IMPROVE YOUR CLIENT'S SKIN WITH THE CONTOUR DERMA II

Using a non-surgical skin treatment approach that combines 2 effective technologies: micro-needling and radio frequency to precisely address specific areas of concern. Through radio frequency and micro-needling, this combined system generates controlled dermal heating, elevating the skin's temperature from 55°C to 65°C, revitalising the skin, triggering the production of new collagen and elastin, simultaneously penetrating the deep dermal and sub-dermal layers. This process significantly enhances skin firmness and tautness from within.

The interchangeable micro-needles, plated in gold and non-insulated, ensure deep penetration. As cell turnover increases, they stimulate the body's natural healing process, prompting the production of collagen and elastin within the dermis. This results in facial and body remodelling and contouring, creating a smoother and more desired appearance while safeguarding the epidermis. Moreover, the power and depth of the treatment can be adjusted to minimise any potential discomfort or downtime.

Contour Derma boasts remarkable results in the realms of anti-aging, face lifting, smoothing, addressing skin laxity, rejuvenating the skin, reducing fine lines and wrinkles, and correcting issues such as scarring and stretch marks.



Precise and Secure Fractional Treatment reaching depths of **3.5mm**, coupled with an extra **2mm** thermal profile.

Minimal dermal thermal impact, making it a secure option for skin type VI, with minimal risk of post-inflammatory hyperpigmentation (PIH).





# TREATMENT SPECS & PRINCIPLES

## LCD TOUCH SCREEN

> USER-FRIENDLY AND EFFORTLESS OPERATION

> SIMPLE STEP-BY-STEP SCREEN GUIDE FOR CHOOSING CORRECT TREATMENT AREAS & DEPTHS



## PROFESSIONAL, SAFE & FAST

> SUITABLE ON ALL SKIN TYPES:  
TREATS FITZPATRICK I - IV SAFELY WITH LOW RISK OF POST  
INFLAMMATORY HYPER-PIGMENTATION

> LITTLE TO NO THERMAL DAMAGE TO DERMIS


> DEEP & SAFE FRACTIONAL TREATMENT PENETRATION  
BETWEEN 0.5MM & 3.5MM

> CONVENIENT TREATMENT WITH HANDLE MOUNTED  
BUTTON SWITCH FOR EASY OPERATION

> EFFORTLESS NEEDLE CARTRIDGE REPLACEMENT PROCESS







### 1. INSULATED NEEDLES

- > Equal therapy to epidermis and dermis layer
- > Adjustable & precise depth control
- > Size 0.5-3.5mm

### 2. STEPPING MOTOR TYPE

- > Needle inserts in to the skin smoothly & effortlessly without dragging

### 3. GOLD PLATED NEEDLES

- > High biocompatibility
- > Suits metal allergy clients

### 4. SAFETY NEEDLE SYSTEM

- > Sterilized disposable needle cartridges
- > Suction combined probe for better skin contact

### 5. USER-FRIENDLY HAND PIECE DESIGN



- > 2 different needle cartridge sizes to suite 2 different treatment areas
- > Hand piece trigger button

# WORKING PRINCIPLE

---

## WHAT IS MICRONEEDLING RF?

This innovative treatment involves fractional skin rejuvenation, which effectively triggers collagen production within the deeper dermal layers. By precisely targeting these underlying skin layers, the fundamental architectural components undergo a natural reorganization, fostering a remarkable anti-aging transformation.



Micro-needling represents a combined and innovative approach to skin rejuvenation, harnessing the body's inherent healing mechanisms to elevate skin aesthetics.

While many individuals gravitate towards micro-needling with radio-frequency for its remarkable anti-aging effects, such as enhancing skin texture and minimising wrinkles, it also holds the power to address trauma and surgery scarring, stretch marks plus burned tissues. Through the generation of collagen and elastin induced by the procedure, a gradual fading of scars and pigmentation occurs, culminating in the revitalisation of your client's skin and complexion.

# CLIENT COMFORT & SAFETY ENSURED

---

## COMFORTABILITY

Our machine uses a vacuum suction to help the needles penetrate the skin. This reduces the amount of force that is needed, which makes the procedure less painful.

## INSULATED NEEDLES

The needles are insulated, which means that they only deliver treatment to the Epidermis layer of the skin. This helps to prevent damage to the Dermis layer, which can cause pain and scarring.

## STEPPING MOTOR TYPE

The stepping motor type needle insertion is smoother and causes less pain than the solenoid type. This is because the stepping motor moves the needles in a controlled way, which reduces the risk of shock.

## GOLD-PLATED NEEDLES (25 PINS - 64 PINS)

The gold-plated needles are durable and have high biocompatibility. This means that they are less likely to cause allergic reactions. Patients with metal allergies can use this treatment without worrying about developing contact dermatitis.

## PRECISE DEPTH CONTROL - 0.5 ~ 3.5MM

The needle depth can be precisely controlled in 0.1mm increments. This allows the treatment to be customised to the individual patient's needs.


## SAFETY NEEDLE SYSTEM

This system helps to ensure that the treatment is safe and effective:

- > Sterilised disposable needle tips
- > A red light that indicates when RF energy is being applied

## ELABORATE NEEDLE THICKNESS - MIN: 0MM < MAX: 0.2MM

The needles are thin and easy to penetrate the skin with minimal resistance. This helps to reduce pain and discomfort during the procedure.





# TREATMENT AREAS

---



Forehead

Brows

Canthal Lines

Cheeks

Upper Lip

Abdomen

Jowls

Jaw Line

Neck

Chest

# TREATMENT TIME & PRICE GUIDANCE


## RADIOFREQUENCY MICRO NEEDLING

TREATMENT AREA	DURATION (mins)	PRICE (£)	COURSE OF 3 (10% saving)	VASCULAR TREATMENTS
<b>Eyebrow Lift</b>	<b>10</b>	<b>50.00</b>	<b>135.00</b>	<b>£150-£500 varies depending on size of treatment area(s).</b>
<b>Forehead</b>	<b>20</b>	<b>70.00</b>	<b>189.00</b>	
<b>Neck</b>	<b>30</b>	<b>70.00</b>	<b>189.00</b>	
<b>Jawline</b>	<b>30</b>	<b>70.00</b>	<b>189.00</b>	
<b>Half Face</b>	<b>45</b>	<b>170.00</b>	<b>459.00</b>	
<b>Full Face</b>	<b>60</b>	<b>250.00</b>	<b>675.00</b>	
<b>Full Face &amp; Neck</b>	<b>75</b>	<b>300.00</b>	<b>810.00</b>	
<b>Body (15 x 15 cm)</b>	<b>45</b>	<b>400.00</b>	<b>1080.00</b>	
<b>Body (30 x 30 cm)</b>	<b>60</b>	<b>700.00</b>	<b>1890.00</b>	

Treatment times are approximate and are given as guidance for the time the machine is in use. (Not full appointment time, as we recommend adding 15 minutes per area for appointment set up, clean down and client interaction.)

## POTENTIAL EARNINGS ARE HIGH!

30 MINUTE TREATMENT <b>£300</b> 60 MINUTE TREATMENT <b>£600</b>	<b>£2,000+</b> PER WEEK	<b>£8,000+</b> PER MONTH	<b>£100,00+</b> PER YEAR
--	----------------------------	-----------------------------	-----------------------------



# OUR PAYMENT OPTIONS



Contour Group offers **easy and affordable purchasing options that perfectly fit your business' cash flow**. We are here to assist you in providing the finest beauty/laser treatments at your clinic or beauty salon, using the latest and most advanced beauty devices available from us.

## INCREASE YOUR PROFITS!

We assist you in unlocking your business' potential, allowing you to boost revenue streams and ultimately generate additional income to fuel the growth of your business! Our cost-effective payment options provide a dual advantage for your clinic, tailored specifically to your unique individual or business circumstances.

### OUR MULTIPLE PAYMENT OPTIONS



#### **PURCHASE OUTRIGHT (SELF-FUNDING)**

If you wish to make a full payment for a machine without any financing options. To initiate the process, [please email us at sales@contourgroup.co.uk](mailto:sales@contourgroup.co.uk) stating the specific machine(s) you are interested in purchasing. We will issue an invoice for payment.



#### **FLEXIBLE FINANCE**

Secure financing through our dedicated partners at Admiral Leasing. Enjoy the flexibility of fixed repayment terms (2-5 years) with or without a deposit. You also have the choice of deferred payments, where you don't need to make any payments for 6 months, or pay just £50.00 for the first 6 months. Select the repayment duration that best aligns with your requirements. (Subject to a finance credit check) **Tel: 0161 628 1250 and ask for our personal reps Carol, Hayley or Frances.**



#### **INTEREST FREE FINANCE**

Pay a required deposit and spread the remaining balance over 6-12 months, on our interest free monthly repayment agreement.



#### **GOVERNMENT BACKED BUSINESS STARTUP LOAN SCHEME**

**An unsecured loan** for new business' or business' trading less than 36 months. Terms available: 1-5 years. Loans available for up to £25k. With no early repayment penalties and a government fixed rate of 6%. [Click here to apply](#)



# ADDITIONAL SERVICES & AFTER SALE SUPPORT

---

We go beyond the purchase of your new machine! Contour Group are proud to provide these additional **FREE** services and benefits to our valued customers:



**BUSINESS LAUNCH MANAGER ALLOCATED FOR 1:1 SUPPORT AND GUIDANCE**



**TREATMENT PRICE GUIDES AND TEMPLATES PROVIDED**



**COMPLETE AFTERSALES SUPPORT**



**MACHINE DELIVERY AT YOUR OWN UK CLINIC LOCATION WITHIN 2-3 WEEKS**



**FULL MACHINE SET UP AND INSTALLATION BY DEDICATED ENGINEER**



**FREE ONLINE ACCREDITED CORE OF KNOWLEDGE LASER TRAINING (VIA PISTACHE-TRAINING.COM)**



**FREE CERTIFICATED MANUFACTURERS MACHINE TRAINING FOR UP TO 5 STAFF MEMBERS IN YOUR OWN UK CLINIC LOCATION**



**ACCESS TO OUR WHATSAPP BUSINESS SUPPORT & HELPFUL HINTS/TIPS GROUP**



**ASSISTANCE WITH MARKETING STRATEGIES TO ADVERTISE YOUR NEW MACHINE TREATMENT**



**ACCESS TO OUR CANVA FILE FOR ADVERTS, PROMOTIONAL MATERIALS & SOCIAL MEDIA GRAPHICS**



**TECHNICAL SUPPORT: UK EXPERIENCED LASER & BEAUTY MACHINE ENGINEERS FOR TROUBLE SHOOTING, REPAIRS & ANNUAL SERVICE REQUIREMENTS**




**ENGINEER VIDEOS AVAILABLE FOR SIMPLE TUTORIALS (I.E. CHANGING WATER, CHANGING FILTERS ETC.)**



**REFERRAL TO OUR INDEPENDENT BRANDING SPECIALIST TO ASSIST WITH EXPERT BRANDING/LOGO/CORPORATE IDENTITY/PROMOTIONAL ADVERT CREATIONS/SOCIAL MEDIA ACCOUNT SETUP & MANAGEMENT. (ADDITIONAL FEES APPLICABLE FOR THESE EXTRA SERVICES)**

# RESULTS ARE EXCELLENT!



## WE ARE CONTOUR GROUP™

Providing the latest cutting edge, innovative and advanced treatment machine solutions for beauty and aesthetic clinics.

The origin of the Contour Group™ was as a salon and beauty clinic looking to make more money. When researching ways to increase staff efficiency and boost overall monthly revenue, the group researched and cultivated excellent machinery to make beauty clinics dramatically increase their profits whilst providing unparalleled treatment results for their clients.

All of our machines come with a 12-month parts and labour warranty.

To book a free Zoom video machine demo, or a face-to-face machine demo, please email:- [sales@contourgroup.co.uk](mailto:sales@contourgroup.co.uk)

